

Troubleshooting guide for myopic orthokeratology with OMEGA CL contact lenses

The following indications must be considered exclusively as a supplement to the practitioner's professional knowledge and skills. The issues often have multifactorial causes.

VISION ISSUES

At distance, neutralized with a negative sphere	Wait a few more days Increase "T" parameter value (more negative)
At distance, does not neutralize with a negative sphere	Epithelial defect Higher order aberrations Residual astigmatism: <ul style="list-style-type: none">• Concurrent corneal toricity:<ul style="list-style-type: none">○ Switch to OMEGA TORIC geometry○ Switch to OMEGA DOUBLE TORIC geometry (technical assistance)• Internal astigmatism:<ul style="list-style-type: none">○ Switch to OMEGA DOUBLE TORIC geometry (technical assistance)
At distance, neutralized with a positive sphere	Reduction of parameter "T" value (less negative) Lens adhesion
At distance, worsening vision in the evening	Scotopic conditions: <ul style="list-style-type: none">• Pupillary diameter• Treatment area diameter• Decentered treatment Treatment regression: <ul style="list-style-type: none">• Increase in parameter "T" value (more negative)
Worse upon awakening	Ocular surface problem Edema
At near, neutralized with positive	Overcorrection: <ul style="list-style-type: none">• Reduction of parameter "T" value (less negative) Glasses for near
At near, does not neutralize with positive	Multifactorial etiology

TREATMENT AREA ISSUES

Downward decentralization

Increase the value of parameter "C"
Switch to OMEGA TORIC geometry

If after a long period of wearing:

- Check lens surface quality and lens care regimen
- Avoid prolonged use with open eyes

Upward decentralization

Decrease the value of parameter "C"
Switch to OMEGA TORIC geometry

Increase the total diameter

If after a long period of wearing:

- Check lens surface quality and lens care regimen
- Avoid prolonged use with open eyes

Decentralization towards the nose / temple

Decrease the value of parameter "C"
Switch to OMEGA TORIC geometry

Increase the total diameter

If after a long period of wearing:

- Check lens surface quality and lens care regimen
- Avoid prolonged use with open eyes

Switch to OMEGA ASYMM (technical support)

Treatment area too small (difficult to resolve)

Increase the value of parameter "C"

Decrease height of OZ

Central island in topography

Check cornea epithelial integrity

Increase height of OZ

FITTING ISSUES

Excessive central pooling

Increase the value of parameter "C"
Decrease height of OZ

Excessive central bearing

Reduce the value of parameter "C"
Increase height of OZ

Incomplete alignment

Reduce the value of parameter "C"
Evaluate OMEGA TORIC geometry

Excessive movement

Excess tearing, wait
Reduce the value of parameter "C"
Increase the total diameter
Evaluate OMEGA TORIC geometry

Limited movement

Increase the value of parameter "C"
Reduce the total diameter
If OMEGA TORIC, reduce toricity

Limited fluorescein penetration under the lens

Increase the value of parameter "C"
Reduce the total diameter
If OMEGA TORIC, reduce toricity

Presence of bubbles under the lens

Improve application proficiency
Apply lens with more viscous solution
If alignment is incomplete, consider OMEGA TORIC geometry

Difficulty with removal

Check for signs of lens adhesion
Increase the value of parameter "C"
Reduce the diameter
If OMEGA TORIC, reduce toricity

OCULAR SURFACE ISSUES

Central epithelial staining	<p>Excessive central bearing:</p> <ul style="list-style-type: none">• Reduce the value of parameter "C"• Increase height of OZ <p>Deposits on the lens:</p> <ul style="list-style-type: none">• Check lens care regimen <p>Epithelial fragility:</p> <ul style="list-style-type: none">• Wait a few nights• Apply lens with more viscous solution
Diffused punctate epithelial staining	<p>Check maintenance procedures / products</p>
Epithelial staining in the area of alignment	<p>Check lens care regimen</p> <p>Increase the value of parameter "C"</p> <p>Reduce the total diameter</p> <p>If OMEGA TORIC, reduce toricity</p>
Impression of lens	<p>Increase the value of parameter "C"</p> <p>Reduce the total diameter</p> <p>If OMEGA TORIC, reduce toricity</p>
Corneal infiltrates	<p>Discontinue treatment and manage condition</p>
Sharp red eye (particularly if monocular)	<p>Discontinue treatment and manage condition</p>
Occasional discomfort	<p>Foreign body</p> <p>Damaged lens</p> <p>Check lens surface quality and lens care regimen</p>